



Transformation of Chemically Competent Cells

1. Chill either 14 ml polypropylene tube or 1.5 ml microfuge tube on ice.
2. Thaw competent cells on ice.
3. Warm SOC to 37°C.
4. Aliquot 50 ul of competent cell per tube
5. Add 1 ul to 5 ul DNA (as appropriate).
6. Incubate on ice 15 minutes.
7. Heat pulse at 42°C for 30 seconds.
8. Return tube to ice for 1 minute.
9. Add anywhere between 200 ul and 950 ul SOC recovery media to tube.
10. Incubate tube in 37°C shaker (225 rpm) for 60 minutes.
11. Plate appropriate volume on selective plate.